



SCCT2023 | BOSTON

July 27 - 30 | Abstract Submission Site

[Print this Page for Your Records](#)

[Close Window](#)

Control/Tracking Number: 2023-A-458-SCCT

Activity: Abstract

Current Date/Time: 3/16/2023 12:32:19 PM

AI-enabled Cardiac Chambers Volumetry In Non-contrast Coronary Artery Calcium CT Scans Vs. Contrast-enhanced Coronary CT Angiography Scans In The Same Patients

Author Block: Chenyu Zhang, MS¹, Kyle Atlas, BS¹, Song Shou Mao, MD², Anthony P. Reeves, PhD³, Matthew Budoff, MD², Claudia Henschke, MD⁴, David Yankelevitz, MD⁴, Edward Callahan, MD¹, Morteza Naghavi, MD¹.

¹HeartLung Technologies, Torrance, CA, USA, ²The Lundquist Institute, Torrance, CA, USA, ³Cornell University, Ithaca, NY, USA, ⁴Mount Sinai Hospital, New York, NY, USA.


Abstract:

Introduction: We have developed an AI-enabled automated volumetry of cardiac chambers (AutoChamber) that works both on non-contrast ECG-gated cardiac CT scans used for coronary artery calcium (CAC) score, and non-contrast non-gated lung cancer screening CT scans. We have recently reported the agreement between AutoChamber measurements vs cardiac MRI measurements. Here we report AutoChamber results on non-contrast CAC scans and contrast enhanced coronary CT angiography (CTA).

Methods: We have obtained CT images from a previously published study in which 131 cases underwent non-contrast CAC and contrast enhanced CCTA at the Harbor UCLA Lundquist Research Institute. The cases consisted of 74 [52.4%] women with a mean age of 55 years \pm 10.1. The cardiac chambers volume in 131 cases of contrast enhanced CCTA were measured using automated chamber volumetry by Philips Brilliance Volume v 4.5; Philips Healthcare.

Results: The correlations between Philips' automated chamber volumetry in contrast enhanced CCTA scans versus AutoChamber automated volumetry in non-contrast CAC scans for left atrium (LA), left ventricle (LV), right atrium (RA), right ventricle (RV), and left ventricular wall (LVW) volume were 0.73, 0.65, 0.85, 0.93, 0.90, respectively. P for all <0.01.

Conclusions: AI-enabled AutoChamber volumetry based on non-contrast cardiac scans is feasible, and correlates reasonably well with contrast-enhanced cardiac CT scans. Further studies are warranted to reproduce these results and explore the add-on value of AutoChamber to CAC scans.



Category (Complete): LV/RV Function, Chamber Dimensions ; Artificial Intelligence/Machine Learning

Abstract Type (Complete):

* **Would you like to be considered for the YIA program?:** N/A

Status: Complete

Society of Cardiovascular Computed Tomography (SCCT)

Telephone: 703.766.1706

e-mail: education@scct.org

Helpdesk:

Phone: 217-398-1792

Email: scct@support.ctimeetingtech.com

💡 [Feedback](#)

Powered by [cOASIS](#), The Online Abstract Submission and Invitation System SM
© 1996 - 2023 [CTI Meeting Technology](#). All rights reserved. [Privacy Policy](#).